



Two Striper.

Compo-Disc®

HAND-HELD DIAMOND FINISHING INSTRUMENT

Use Compo-Disc® diamond discs to open proximal contacts and for contouring. Diamond Finishing Strips are ideal for smoothing and finishing all proximal restorations and veneer margins.

Compo-Disc® ORDER NUMBER: CECSTFTP106

- Thin
- Edge and face-cutting
- Single-sided
- Patented holder design

Two Striper.

Diamond Finishing Strips

HAND-HELD DIAMOND FINISHING INSTRUMENT

Diamond Finishing Strips

ASSORTED PACK
ORDER NUMBER: CEC50087

- No diamond stripping or peeling
- Color Coded
- Flexible
- Safe-sided
- Sterilizable



ORDER NUMBER	WIDTHS	THINNESS	COLOR	GRIT
100T	2.5 MM (0.10 IN)	0.15 MM (0.006 IN)	BLUE	60μ
150T	3.75 MM (0.15 IN)			
100F	2.5 MM (0.10 IN)	0.127 MM (0.005 IN)	RED	45μ
150F	3.75 MM (0.15 IN)			
100UF	2.5 MM (0.10 IN)	0.1 MM (0.004 IN)	YELLOW	20μ
150UF	3.75 MM (0.15 IN)			

Compo-Disc® - Package of 2 discs Diamond Finishing Strips - Package of 6 or 12 Assorted Strips (Also available in individual packages of 6 or 12 for each strip)

Diamond Finishing Strips 200T2 - Package of 3

Diamond Finishing Strips 200T2

ASSORTED PACK
ORDER NUMBER: CERPOT3P200

For Orthodontic Reproximation

Use the Diamond Finishing Strips 200T2 for interproximal enamel reduction.

- Flexible, double sided
- 60 micron diamond
- Safe center
- Fast and smooth cutting
- Sterilizable
- Ultrasonic compatible